

Advancing High Throughput Screening of CFTR Modulators with Human Epithelial Cell Lines

Jinliang Sui, Haibo Shang, Abhijeet Kanawade, Richard Fitzpatrick
Flatley Discovery Lab, LLC. 529 Main Street, Charlestown, MA 02129

ABSTRACT

Human epithelial cell lines were adapted to the iodide flux HTS assay: A549 cells stably transfected with YFP and CFTR-F508del (Pedemonte et al., 2010, AJP 298), and CFBE410-transiently transfected with YFP and CFTR-F508del. For CFBE410- cells, a suitable stably transfected cell line has not been available. We introduced a high volume transient transfection method (STX by MaxCyte) for the CFBE410- cells to express YFP and CFTR-F508del for the HTS. This method allows plasmid titration in order to achieve desirable expression levels and expression ratios of YFP and CFTR-F508del for best assay signals, with good transfection efficiency and cell viability. This method also demonstrated consistency across experiments and for various freezing durations. The high transfection volume and the consistent transfection quality make this method a suitable approach for the HTS. Data generated from FRT, A549, and from CFBE cells against a set of selected compounds are compared, showing the benefits of using human epithelial cells in the HTS assay. Advancing the iodide flux assay on the HTS platform from Fisher Rat Thyroid cells to human epithelial cells is a significant step towards the ideal CF cellular model. The human cell-based iodide flux assays are currently the primary HTS methods in the CFTR modulator screening campaign at Flatley Discovery Lab.




Figure 1. Sample microscope images of transfected CFBE, A549, and FRT cells in phase-contrast view (top panel), GFP view (middle panel), and merged view (lower panel), showing the expression of YFP in the cells under the conditions where they are ready for the experiments




Figure 2. Sample records of YFP signals from FRT cells (A), CFBE cells (B), and A549 (C) cells, showing traces from negative control wells (cells were treated with DMSO only, DMSO) and from positive control wells (cells were treated with C18 for 24-hours, C18) respectively.




Figure 3. Comparison of YFP quenching rate (QR) signal windows between negative control (DMSO) and positive control (C18, 10 um) among the three tested cell lines (FRT, A549, and CFBE) in the YFP assay





Figure 6. CFTR Corrector validation hits

Table data is sorted by ADX of CFBE cell data, from a sample data set from the iodide flux HTS assay with three different cell lines (FRT, A549, and CFBE410-), showing the top performers in the assay including VX-809, C18, C17, SAHA, C4.




Figure 5. CFBE410- iodide flux assay performance in one experiment and across multiple experiments.

A. plot of differential quenching rate (dQR) and z-factors from one experiment set of assay plates with CFBE410- cells expressing YFP & CFTR-F508del transiently transfected with MaxCyte STX; B. z-factor plot of 672 assay plates across multiple experiments.

MaxCyte® STX™

Scalable Transfection System

Proprietary Electroporation Technology

Simple:

High Yield:

High Efficiency:

Safe:

Scalable & Rapid:

Rugged:

Quality:

Fast and easy

>90% viability

>90% transfection efficiency

Chemically defined buffer

No added biological agents

Sterile, closed system

5 x 10⁵-4x10⁷ (Static EP) in seconds

1 x 10¹⁰ cells (Flow EP) in <30 min

Reproducible & consistent

cGMP compliant

ISO 9001 certification

CE Marking

Small molecules

Antigens (proteins/lysates)

Nucleic acids (DNA, mRNA, siRNA)

Primary cells

Stem cells

Mammalian cell lines

Figure 7. The MaxCyte STX system (courtesy of James Brady from MaxCyte Inc.)

Figure 8. Large Scale transfection of GFP with CFBE410- cells EP protocol = "A549". [DNA] = 100 ug/mL, [cell] = 3e7/mL. FACS Data 24 hrs post electroporation (courtesy of James Brady from MaxCyte Inc.).

INTRODUCTION

The YFP-based iodide influx assay has been widely deployed in high-throughput screening (HTS) for small molecule CFTR modulators. It's a robust cell-based assay, developed with Fisher Rat Thyroid (FRT) cells expressing YFP and CFTR-F508del (Pedemonte et al., 2005, JCI 115; Sui et al., 2010, ADDT 8).

It has been highly desirable to use human epithelial cells replacing FRT cells in this assay to eliminate concerns over the species gap between rat and human. It may significantly boost our chance of finding clinically relevant HTS hits in our efforts seeking small molecule CFTR modulators. This work is to report our work along this line using A549 and CFBE cells, both of which are human bronchial cell lines, to replace FRT cells.

MATERIALS & METHODS

FRT cells and A549 cells stably expressing YFP and CFTR-F508del were kindly provided by Dr. Galiotta. High volume transient transfection for HTS was developed with MaxCyte STX (see pictures 7 & 8). CFBE410- cells were transfected and frozen immediately. Freshly thawed cells were seeded in assay plates and treated with testing compounds following the CFTR corrector protocol for the iodide influx assay (Sui et al, ADDT 2010).

YFP data analysis, HTS data quality control, and the HTS and chemical library data management software system were developed in-house. YFP signal quenching rate (QR) was derived by a non-linear curve fitting to the equation $f(t) = (A_0 - A_1)e^{-tQR} + A_1$

Differential quenching rate (dQR) is given by $dQR = QR - QR_{NC}$, where QR_{NC} is the mean QR of the DMSO control wells in the assay plate. For each compound with the results of all three (0.2μM, 2μM, and 20μM) or five (0.2μM, 0.66μM, 2μM, 6.6μM, and 20μM) testing doses, the activity index (ADX) is calculated with the formula $ADX = \sum_{k=1}^n (A_k - A_{k-1}) [-\log(C_k)]$, where A_k is the activity measurement, C_k is the testing concentration.

CONCLUSION

• At Flatley Discovery Lab (FDL), we have adopted several human epithelial cell lines for the iodide flux HTS assays.

• A high volume transient transfection method (STX by MaxCyte) has been established to express YFP and CFTR-F508del in the human bronchial cells for the HTS assays. This method allows plasmid titration in order to achieve desirable expression levels, with consistently high transfection efficiency and cell viability.

• In the CFTR modulator screening campaign at FDL, we are using the human bronchial cells in the iodide flux HTS primary assays, and we have screened several hundreds of thousands of library compounds. This effort has yielded several interesting hits.

Acknowledgements

We thank Dr. Luis Galiotta at Gaslini Institute, Genova, Italy for his kind contribution of the A549 cells; We are also grateful to Dr. Martin Mense and Herman Bihler at CFFT for their support and participation in this work. This work is supported by the Flatley Foundation.

References

Pedemonte N, Tomati V, Sondo E, Galiotta LJ. Influence of cell background on pharmacological rescue of mutant CFTR. Am J Physiol Cell Physiol. 2010 Apr;298(4):C866-74.
Sui J, Cotard S, Andersen J, Zhu P, Staunton J, Lee M, Lin S. Optimization of a Yellow fluorescent protein-based iodide influx high-throughput screening assay for cystic fibrosis transmembrane conductance regulator (CFTR) modulators. Assay Drug Dev Technol. 2010 Dec;8(6):656-68. Epub 2010 Nov 4.

Flatley Discovery Lab

Seeking A Cure for Cystic Fibrosis